

Delivery Devices

BONASTENT®

- * Superior Pushability and Flexibility
- * Kinking Resistance
- * Design for Easy & Accurate Placement

Y - Shape Handle with Braided Tubing

10Fr or Less Delivery Device for

- Biliary Non-Covered & Covered
- Pyloric / Duodenal Non-Covered & Covered
- Colorectal Non-Covered & Covered
- Tracheal / Bronchial



Radiopaque Band

Radiopaque Marking for Distal End of Delivery Device



Yellow Markin Band

Endoscopic View Marking for Proximal End of Captured Stent



I - Shape Handle with Teflon Tubing

Larger than 10Fr Delivery Device for

- Esophageal & Esophageal with Shim Anti-Reflux Valve
- Colorectal Non-Covered & Covered
- Tracheal / Bronchial



Centimeter Marking

Positioning Measurement for Esophageal & A.R.V. Stent